



### Quick Preset

With one touch, the user can select the most common transducer and preset combinations. Quick Preset maximizes efficiency to make a full day of scanning simple and easy.




### Touch Customization

Samsung has made a customizable touchscreen interface that allows the user to move frequently used functions to the first page, keeping the focus on the patient instead of the system.



### 6-way Control Panel

The RS85's 6-way adjustable control panel optimizes your work environment to reduce repetitive motions stress. When it's in off-mode, the control panel returns to the home position, allowing for easier and enhanced mobility.



### 13.3-inch Tilting Touch Screen

Samsung's tilting touch screen can be adjusted to accommodate any user's viewing preferences within any scanning environment.



### Central Lock

A single pedal controls a central lock mechanism to conveniently secure the console in place. This results in more efficient movements while the user is performing scanning procedures.



### Gel Warmer

Samsung's two-level adjustable gel warmer keeps ultrasound gel at a comfortable temperature.



### Maneuverable Wheel

4 swivel wheels allow easy steering, and a locking function.

