



QuickPreset

With one touch, the user can select the most common transducer and preset combinations. QuickPreset increases efficiency to make a full day of scanning simple and easy.




EzCompare™

EzCompare™ allows easy access to previously taken exams to evaluate corresponding views in a side-by-side display. For greater efficiency, EzCompare™ automatically matches the image settings, annotations, and bodymarkers from the prior study.



Measure Navigation

When placing a caliper, Measure Navigation automatically magnifies the region of interest using a picture-in-picture window to allow more precise placement of the calipers. This is especially useful when measuring small structures or when accuracy is critical.





Gel warmer

For operator convenience, a gel warmer can be installed on both sides of the control panel.



Solid State Drive (SSD)

The HS50 uses Samsung's advanced solid state drives. These stable and dependable drives allow faster boot-up, better frame rates, and fast processing speeds.



BatteryAssist™

BatteryAssist™ provides the system with battery power. This serves two important purposes. Firstly, it enables users to perform scans and transport the ultrasound system to other locations in environments where AC power may not be available temporarily. Secondly, it reduces boot-up time by using sleep mode without having to shut down or restart the system.



Clever use of space

With its reduced weight and compact size, the HS50 takes up minimal space and can move freely. In addition, its streamlined rear profile allows you to park the HS50 in small spaces.



Rear Tray

HS30's rear tray provides extra storage space for the endo-cavity transducer and other items.



Gel Warmer

Two-level adjustable gel warmer maintains ultrasound gel at a comfortable temperature.



Side Pocket

The system side pocket is ideal for storing a tablet, patient charts, or other items that you need to keep close at hand.



21.5-inch LED Monitor

21.5-inch full HD LED display delivers excellent contrast resolution, image clarity and vibrant color.



Keyboard & Keyskin

The soft typing designed keyboard helps quick actuation and optimizes working environment. Keyskin for the keyboard and control panel is helpful to avoid any contaminants.



Transducer Cable Hangers

Users can arrange the transducer cables neatly on the 2 hangers on either side of the system.